









ONE-CARE

GLUCOSE MONITORING SYSTEM
FOR SELF-TESTING

Glucose monitoring has never been easier...



 FREE Auto Calibration	 0.5µl 0.5 µL Blood Sample	 AST Alternative Site Testing
 5 sec 5 Second Test	 500 500 Memory Sets	 PC Download
 14 Instant 14 Days Average	 Strip Ejector	



Bigger Buttons in a compact device

Easy to use for elderly people



USB and Bluetooth (Optional)

Provide both wire and wireless connection for data transmission



Test Strip Ejector

Reduce the possibility of cross-contamination due to bloodborne pathogens



Test Strip

Longer shelf life 24 months



ONE-CARE

SPECIFICATIONS




Auto Coding

Simplified way to test, easy to use, no more misreading due to incompatible code card

GOD Enzyme

0.5 µl blood sample
Highly Accurate for Capillary Whole Blood samples




ONE-CARE™

- Hematocrit** | 30-55%
- Test Sample** | Capillary Whole Blood
- Measuring Unit** | mg/dL or mmol/L
- Measuring Range** | 20 - 600 mg/dL (1,1 - 33,3 mmol/L)
- Average** | 14/07/30 days
- Reminder Alarm** | 4 sets
- External Output** | UART with mini USB sock (Optional)
| Bluetooth (Optional)
- Power Supply** | 2 AAA batteries
- Battery Life** | Approximately 1,000 tests
- Dimensions** | 100 (L) x 57 (L) x 24 (A) mm
- Weight** | 58 g without batteries


INSTRUCTIONS OF USE

1 STEP




Insert a Test Strip

2 STEP




Apply Blood Sample

3 STEP



Wait 5 Seconds
and Check Results

4 STEP



Discard the used
Test Strip

ORDERING INFORMATION

TEST KIT



TEST KIT

- > 1 Blood Glucose Meter with batteries
- > 10 Blood Glucose Test Strips
- > 1 Adjustable Depth Lancing Device
- > 10 Sterilized Lancets
- > 1 Protective Carrying Case
- > User Manual, Package Inserts

Note: The contents of a Complete Kit may vary and are subject to change according to the actual purchase.

CODE	PRODUCT	PCS per BOX	BOXES per CASE
4006	ONE-CARE® GLUCOSE MONITORING KIT FOR SELF-TESTING	1	80
4007	ONE-CARE® GLUCOSE STRIPS FOR SELF-TESTING	50	560
4010	ONE-CARE® Control Solution Level 1 & 3 (Low & High)	2 bottles	900 boxes (1800 bottles)

WWW.ONE-CARE.COM

MediVena®
inspired by innovation

Americas

USA
MediVena
1 Glenlake Parkway, Ste 700
Atlanta, GA 30328, USA
+1 404 422 1513
info@medivena.com
www.medivena.com

Europe

Poland
MediVena Sp. z o.o.
Al. Jana Pawła II 27
00-867 Warsaw, Poland
+48 695 995 700
info@medivena.com
www.medivena.com

