

# ONE-CARE®

GLUCOSE MONITORING SYSTEM  
FOR SELF-TESTING

Glucose monitoring has never been easier...



Auto Calibration	0.5 µL Blood Sample	Alternative Site Testing	Instant 14 Days Average
5 Second Test	500 Memory Sets	PC Download	Strip Ejector



## Bigger Buttons in a compact device

Easy to use for elderly people



## USB and Bluetooth (Optional)

Provide both wire and wireless connection for data transmission



## Test Strip Ejector

Reduce the possibility of cross-contamination due to bloodborne pathogens



## Test Strip

Longer shelf life 24 months



SPECIFICATIONS




**Auto Coding**

Simplified way to test, easy to use, no more misreading due to incompatible code card

**GOD Enzyme**


0.5 µl blood sample  
Highly Accurate for Capillary Whole Blood samples



Hematocrit	30-55%
Test Sample	Capillary Whole Blood
Measuring Unit	mg/dL or mmol/L
Measuring Range	20 - 600 mg/dL (1,1 - 33,3 mmol/L)
Average	14/07/30 days
Reminder Alarm	4 sets
External Output	UART with mini USB socket (Optional)   Bluetooth (Optional)
Power Supply	2 AAA batteries
Battery Life	Approximately 1,000 tests
Dimensions	100 (L) x 57 (L) x 24 (A) mm
Weight	58 g without batteries


INSTRUCTIONS OF USE

**1 STEP**




Insert a Test Strip

**2 STEP**




Apply Blood Sample

**3 STEP**



Wait 5 Seconds and Check Results

**4 STEP**



Discard the used Test Strip

ORDERING INFORMATION

TEST KIT



TEST KIT

- > 1 Blood Glucose Meter with batteries
- > 10 Blood Glucose Test Strips
- > 1 Adjustable Depth Lancing Device
- > 10 Sterilized Lancets
- > 1 Protective Carrying Case
- > User Manual, Package Inserts

**Note:** The contents of a Complete Kit may vary and are subject to change according to the actual purchase.

CODE	PRODUCT	PCS per BOX	BOXES per CASE
4006	ONE-CARE® PRO GLUCOSE MONITORING KIT FOR SELF-TESTING	1	80
4007	ONE-CARE® PRO GLUCOSE STRIPS FOR SELF-TESTING	50	560

WWW.ONE-CARE.COM

Americas

**USA**  
MediVena  
1 Glenlake Parkway, Ste 700  
Atlanta, GA 30328, USA  
+1 404 422 1513  
info@medivena.com  
www.medivena.com

Europe

**Poland**  
MediVena Sp. z o.o.  
Al. Jana Pawła II 27  
00-867 Warsaw, Poland  
+48 22 299 44 19  
info@medivena.com  
www.medivena.com

