

ONE-CARE[®]

INSULIN & TB SAFETY SYRINGES

OPTI



Your Ideal Safety Syringes with Fixed Needles!


- **Convenient**

Easy pull-out

- **Comfortable**

Sharp and lubricated needles for gentle insertion

- **One-Step and One-Hand Activation**



- **LDV (Low Dead Volume)**

Designed to reduce medical waste (<0.035ml)

- **Precise graduation**


marking for improved visibility and precise dosing

- **Thin-wall syringe barrel**

for smooth injection and reduced consumption of plastic

INSTRUCTIONS FOR USE

- Intuitive handling with simple steps.



ORDERING INFORMATION

- Available in different syringe sizes:

	CODE	SIZE	NEEDLE (DIAMETER x LENGTH)	PACKAGING
INSULIN SAFETY SYRINGES OPTI	2141	1 mL	29G x 1/2" 0.33mm x 12.7mm	100 pcs / bx 8 bx (800 pcs) / cs
	2142	0.5 mL	29G x 1/2" 0.33mm x 12.7mm	100 pcs / bx 8 bx (800 pcs) / cs
TUBERCULIN SAFETY SYRINGES OPTI	2143	1 mL	28G x 1/2" 0.36mm x 12.7mm	100 pcs / bx 8 bx (800 pcs) / cs

WWW.ONE-CARE.COM

MediVena[®]
inspired by innovation

www.medivena.com

AMERICAS:
MediVena
1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA
+1 404 514 2586 info@medivena.com
CUSTOMER SERVICE TOLL FREE: 888 270 6457

EUROPE:
MediVena Sp. z o.o.
Al. Jana Pawła II 27, 00-867 Warsaw, POLAND
+48 695 995 700

