

ONE-CARE[®]

SAFETY WINGED BLOOD COLLECTION SETS

ERGO-SHIELD



Ergonomic and cost-effective solution ...

DEPF-free, low memory and flexible tubing

assists in handling and preserves the blood chemical and physiological composition

Transparent Tube Holder

for better visibility of collection tubes filling up with blood, tube holders provided separately or pre-attached



Tri-bevel, ultra-sharp, electropolished and lubricated needle

for patient comfort



Ergonomic wing design

permits improved positioning during application



Ergonomic Safety Shield

enables one-hand activation of safety mechanism, audible click confirms that needle is securely locked









ONE-CARE[®]

SAFETY WINGED BLOOD COLLECTION SETS
LUER ADAPTER **ERGO-SHIELD**

INSTRUCTIONS FOR USE

- Intuitive handling with simple steps.



ORDERING INFORMATION

- Available in different needle configurations and with/without pre-attached tube holders:

Without Holders

REF	PRODUCT	COLOR	PACKAGING
5050	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 25G x 3/4" , Luer Adapter, 12" Tube	Orange	100 pcs / box 800 pcs / case
5051	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 23G x 3/4" , Luer Adapter, 12" Tube	Blue	100 pcs / box 800 pcs / case
5052	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 21G x 3/4" , Luer Adapter, 12" Tube	Green	100 pcs / box 800 pcs / case

With Tube Holders

REF	PRODUCT	COLOR	PACKAGING
5050-H	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 25G x 3/4" , Pre-Attached Holder, 12" Tube	Orange	50 pcs / box 400 pcs / case
5051-H	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 23G x 3/4" , Pre-Attached Holder, 12" Tube	Blue	50 pcs / box 400 pcs / case
5052-H	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 21G x 3/4" , Pre-Attached Holder, 12" Tube	Green	50 pcs / box 400 pcs / case

WWW.ONE-CARE.COM

MediVena®
inspired by innovation

www.medivena.com

AMERICAS:
MediVena
1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA
+1 404 514 2586 info@medivena.com
CUSTOMER SERVICE TOLL FREE: 888 270 6457

EUROPE:
MediVena Sp. z o.o.
Al. Jana Pawla II 27, 00-867 Warsaw, POLAND
+48 695 995 700

