



ONE-CARE

Blood Collection Systems



Safety Blood Collection Needles
Safety Blood Collection Wing Sets – Glide & Lock
Safety Blood Collection Wing Sets – Press & Lock
Blood Collection Needles
Blood Collection Wing Sets
Tube Holders

ONE-CARE® SAFETY BLOOD COLLECTION NEEDLES



FEATURES

- Size: 18G, 20G, 21G, 22G, 23G
- Injection & Puncture Instrument
- Vacuum Needle
- Instrument Classification Class II

MATERIALS

- ✓ ABS
- ✓ Rubber
- ✓ SUS304 Cannula
- ✓ Silicone Oil

REGULATORY APPROVALS

- CE
- ISO

STERILIZATION / EXPIRY


- ETO, High Temperature Sterilization
- Shelf life: 5 years from mfg date



ORDERING INFORMATION

CODE	PRODUCT	COLOR	PACKAGING
5020	ONE-CARE® Safety Blood Collection Needles – Eclipse without Holder 22Gx1 1/4"	Black	48 / bx 50 bx (2400 ea) / cs
5021	ONE-CARE® Safety Blood Collection Needles – Eclipse without Holder 21Gx1 1/4"	Green	
5022	ONE-CARE® Safety Blood Collection Needles – Eclipse without Holder 19G x 1 1/4"	Cream	
5025	ONE-CARE® Safety Blood Collection Needles – Eclipse with Holder 22Gx1 1/4"	Black	100 / bx 5 bx (500 ea) / cs
5026	ONE-CARE® Safety Blood Collection Needles – Eclipse with Holder 21Gx1 1/4"	Green	
5027	ONE-CARE® Safety Blood Collection Needles – Eclipse with Holder 19G x 1 1/4"	Cream	

ONE-CARE® SAFETY BLOOD COLLECTION WING SETS
GLIDE & LOCK



FEATURES

- Reliable Safety Mechanism for Permanent Lock
- Ergonomic Wings for Ease of Use
- Ultra-sharp Needle for Patient Comfort

ORDERING INFORMATION

	CODE	PRODUCT	COLOR	PACKAGING
Safety Blood Collection Wing Sets Glide & Lock without Holders	5040	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 25Gx3/4"	Orange	50 / bx 20 bx (1000 ea) / cs
	5041	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 24Gx3/4"	Purple	
	5042	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 23Gx3/4"	Blue	
	5043	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 22Gx3/4"	Black	
	5044	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 21Gx3/4"	Green	
	5045	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 20Gx3/4"	Yellow	
	5046	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock without Holders 19Gx3/4"	Cream	
Safety Blood Collection Wing Sets Glide & Lock with Holders	5040-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 25Gx3/4"	Orange	25 / bx 20 bx (5000 ea) / cs
	5041-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 24Gx3/4"	Purple	
	5042-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 23Gx3/4"	Blue	
	5043-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 22Gx3/4"	Black	
	5044-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 21Gx3/4"	Green	
	5045-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 20Gx3/4"	Yellow	
	5046-T	ONE-CARE® Safety Blood Collection Wing Sets – Glide & Lock with Holders 19Gx3/4"	Cream	

ONE-CARE® SAFETY BLOOD COLLECTION WING SETS

PRESS & LOCK



FEATURES & BENEFITS

- Safety assured, needle stick prevention
- Sliding cartridge design, simple and safe
- Ultra-sharp needle minimizes patient discomfort
- More comfortable double wing design, easy operation
- Color coded for easy identification of needle size
- Sterile pack, single use only

REGULATORY APPROVALS

- MDD 93/42/EEC
- US FDA 510K

STERILIZATION / EXPIRY

- EO Sterilized
- Shelf life: 5 years from mfg date

MATERIALS

- ✓ Needle sleeve Polypropylene
 - ✓ Blood collection needle SUS304
 - ✓ Double wing PVC
 - ✓ Safety protective shield MABS
 - ✓ Safety needle handle MABS
 - ✓ Lock connection MABS
 - ✓ Tubing PVC
 - ✓ Connection set ABS
 - ✓ Puncture needle SUS304
 - ✓ Puncture needle hub MABS
 - ✓ Rubber sleeve Isoprene rubber
 - ✓ Puncture needle sleeve Polypropylene
 - ✓ Holder (optional) Polypropylene
- All Materials are latex & DEHP free

ORDERING INFORMATION

ONE-CARE® SAFETY BLOOD COLLECTION WING SETS			
REF	PRODUCT	COLOR	PACKAGING
5060	ONE-CARE® Safety Blood Collection Wing Sets, Press & Lock, 25G x 3/4" Needle, 12"Tube	Orange	50 / bx 4 bx / cs
5061	ONE-CARE® Safety Blood Collection Wing Sets, Press & Lock, 23G x 3/4" Needle, 12"Tube	Blue	
5062	ONE-CARE® Safety Blood Collection Wing Sets, Press & Lock, 21G x 3/4" Needle, 12"Tube	Green	

ONE-CARE® SAFETY BLOOD COLLECTION WING SETS WITH HOLDERS			
CODE	PRODUCT	COLOR	PACKAGING
5060-T	ONE-CARE® Safety Blood Collection Wing Sets, Press & Lock, 25G x 3/4" Needle, 12"Tube with Holder	Orange	50 / bx 4 bx / cs
5061-T	ONE-CARE® Safety Blood Collection Wing Sets, Press & Lock, 23G x 3/4" Needle, 12"Tube with Holder	Blue	
5062-T	ONE-CARE® Safety Blood Collection Wing Sets, Press & Lock, 21G x 3/4" Needle, 12"Tube with Holder	Green	

ONE-CARE® BLOOD COLLECTION NEEDLES



FEATURES

- Screw-on needle hub compatible with all available holders
- High quality lancet point, silicone treated needle tip for smooth penetration
- Thin wall tube to ensure optimal flow
- Sterile, Non-pyrogenic, non-toxic
- Single use only
- Color coding for easier identification of needle size diameter
- 5 year shelf life, Compliant with ISO 7864:2016

MATERIALS



- ✓ Needle hub - Polypropylene + color masterbatch
- ✓ Cannula – Stainless steel 1.4301 (medical grade)
- ✓ Silicone (medical grade)
- ✓ Glue (medical grade)
- ✓ Rubber sleeve – synthetic rubber (LATEX free)
- ✓ Long protective cap – Polypropylene
- ✓ Short protective cap – Polypropylene

ORDERING INFORMATION

REF	PRODUCT	COLOR	PACKAGING
5001	ONE-CARE® Blood Collection Needles 23G x 1"	Blue	100 / bx 10 bx (1000 ea) / cs
5002	ONE-CARE® Blood Collection Needles 23G x 1 1/2"	Blue	
5003	ONE-CARE® Blood Collection Needles 22G x 1"	Black	
5004	ONE-CARE® Blood Collection Needles 22G x 1 1/2"	Black	
5005	ONE-CARE® Blood Collection Needles 21G x 1"	Green	
5006	ONE-CARE® Blood Collection Needles 21G x 1 1/2"	Green	
5007	ONE-CARE® Blood Collection Needles 20G x 1"	Yellow	
5008	ONE-CARE® Blood Collection Needles 20G x 1 1/2"	Yellow	
5009	ONE-CARE® Blood Collection Needles 18G x 1"	Pink	
5010	ONE-CARE® Blood Collection Needles 18G x 1 1/2"	Pink	

ONE-CARE® BLOOD COLLECTION WING SETS



FEATURES

- Size: 0.9x19
- Blood collection needle with holder
- Vacuum Needle
- Instrument Classification Class II

MATERIALS

- ✓ Medical grade PP
- ✓ SUS304
- ✓ ABS

REGULATORY APPROVALS

- CE
- ISO

STERILIZATION / EXPIRY

- ETO
- Shelf life: 5 years from mfg date

ORDERING INFORMATION

	CODE	PRODUCT	COLOR	PACKAGING
Blood Collection Wing Sets without Holders	5030	ONE-CARE® Blood Collection Wing Sets without Holders 25G	Orange	100ea / bx 20 bx (2000 ea) / cs
	5031	ONE-CARE® Blood Collection Wing Sets without Holders 24G	Purple	
	5032	ONE-CARE® Blood Collection Wing Sets without Holders 23G	Blue	
	5033	ONE-CARE® Blood Collection Wing Sets without Holders 22G	Black	
	5034	ONE-CARE® Blood Collection Wing Sets without Holders 21G	Green	
	5035	ONE-CARE® Blood Collection Wing Sets without Holders 20G	Yellow	
	5036	ONE-CARE® Blood Collection Wing Sets without Holders 19G	Cream	
Blood Collection Wing Sets with Holders	5030-T	ONE-CARE® Blood Collection Wing Sets with Holders 25G	Orange	100ea / bx 20 bx (2000 ea) / cs
	5031-T	ONE-CARE® Blood Collection Wing Sets with Holders 24G	Purple	
	5032-T	ONE-CARE® Blood Collection Wing Sets with Holders 23G	Blue	
	5033-T	ONE-CARE® Blood Collection Wing Sets with Holders 22G	Black	
	5034-T	ONE-CARE® Blood Collection Wing Sets with Holders 21G	Green	
	5035-T	ONE-CARE® Blood Collection Wing Sets with Holders 20G	Yellow	
	5036-T	ONE-CARE® Blood Collection Wing Sets with Holders 19G	Cream	

ONE-CARE® TUBE HOLDERS



• Tube holder



• Blood Transfer Device with Female Luer-Lock



• Blood Transfer Device with Male Luer Adapter

FEATURES

- Size: 0.9x19
- Blood collection needle with holder
- Vacuum Needle
- Instrument Classification Class II

MATERIALS

- ✓ Medical grade PP
- ✓ SUS304
- ✓ ABS

REGULATORY APPROVALS

- CE
- ISO

STERILIZATION / EXPIRY

- ETO
- Shelf life: 5 years from mfg date



ORDERING INFORMATION

CODE	PRODUCT	COLOR	PACKAGING
5015	ONE-CARE® Tube Holder	Transparent	100 / bx 20 bx (2000 ea) / cs
5016	ONE-CARE® Safety Tube Holder	Transparent	100 / bx 10 bx (1000 ea) / cs
5017	ONE-CARE® Blood Transfer Device with Female Luer-Lock and Back-End Needle	Transparent	200 / bx 4 bx (800 ea) / cs
5018	ONE-CARE® Blood Transfer Device with Male Luer Adapter and Back-End Needle	Transparent	200 / bx 4 bx (800 ea) / cs

MediVena®

inspired by innovation

Americas

USA

📍 MediVena
1 Glenlake Parkway, Ste 700
Atlanta, GA 30328, USA
☎ +1 404 422 1513
✉ info@medivena.com
🌐 www.medivena.com

Europe

Poland

📍 MediVena Sp. z o.o.
Al. Jana Pawła II 27
00-867 Warsaw, Poland
☎ +48 661 920 142
✉ info@medivena.com
🌐 www.medivena.com

