

Product Catalog 2024

ONE-CARE®



- **Blood Microsampling**
- **I.V. Blood Collection**
- **Glucose Monitoring**
- **Hypodermic Injection Devices**
- **Personal Protective Equipment**

ONE-CARE

SAFETY LANCETS



Ergonomic

Ergonomic body for comfortable grip

Easy

Wide cap for easy twist off



Comfortable















Dual spring safety activation mechanism for enhanced patient's comfort



Every product detail designed with ease of use and comfort in mind

PRODUCT ORDERING INFORMATION

Available in seven color-coded versions to meet different test blood flow requirements and address a variety of patient skin types.

REF	VERSION	DIAMETER (G) / DEPTH (mm/in)	COLOR	PACKAGING
1101	Micro Flow	30G / 1.5 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1102	Mini Flow	28G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1103	Low Flow	26G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1104	Normal Flow	23G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1105	Universal Flow	21G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1106	High Flow	21G / 2.4 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1107	Max Flow <small>Blade</small>	18G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs



For more information direct to your smartphone

SCAN THE CODE
> SAFETY LANCETS





















ONE-CARE[®]

SAFETY LANCETS **OPTI**



PRODUCT SELECTION CHART

Available in ten color-coded versions to meet different test blood flow requirements and address a variety of patient skin types.

REF	DIAMETER (G) / DEPTH (mm/in)	COLOR	PACKAGING
1031	30G / 1.5 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1032	28G / 1.6 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1033	25G / 1.6 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1034	23G / 1.8 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1035	21G / 1.8 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1036	21G / 2.4 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1037	21G / 2.8 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1038	17G Blade / 1.0 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1039	17G Blade / 1.5 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs
1010	17G Blade / 2.0 mm	 	200 pcs / bx 18 bx (3,600 pcs) / cs

ONE-CARE[®]

ADJUSTABLE SAFETY LANCETS **OPTI+**



**AVAILABLE
1Q-2024**

FEATURES

- 3 adjustable depth settings/: 1.4 (Low) | 1.8 (Medium) | 2.2 mm (High)
- Ergonomic & intuitive T-shape design
- Reliable, dual-spring safety mechanism
- 3-facet, siliconized needle for gentle lancing



ORDERING INFORMATION

REF	DESCRIPTION	COLOR	PACKAGING
1050	ONE-CARE OPTI+ Adjustable Safety Lancets Top Button Activated, 23G	Purple	200 pcs / bx 18 bx (3,600 pcs) / cs
1051	ONE-CARE OPTI+ Adjustable Safety Lancets Top Button Activated, 26G	Yellow	200 pcs / bx 18 bx (3,600 pcs) / cs
1052	ONE-CARE OPTI+ Adjustable Safety Lancets Top Button Activated, 28G	Blue	200 pcs / bx 18 bx (3,600 pcs) / cs

ONE-CARE[®]

PLUS

CONTACT-ACTIVATED SAFETY LANCETS

Ergonomic
Ergonomic body for comfortable grip













Reliable
Reliable safety activation mechanism

Comfortable
Wide cap for easy twist off



PRODUCT ORDERING INFORMATION

Available in six color-coded versions to meet different test blood flow requirements and address a variety of patient skin types.

REF	VERSION	DIAMETER (G) / DEPTH (mm/in)	COLOR	PACKAGING
1001	Micro Flow	30G / 1.5 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1002	Mini Flow	28G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1003	Low Flow	25G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1004	Normal Flow	21G / 1.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1005	High Flow	21G / 2.4 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs
1006	Max Flow	21G / 2.8 mm 		100 pcs / bx 20 bx (2,000 pcs) / cs



For more information direct to your smartphone

SCAN THE CODE
➤ PLUS SAFETY LANCETS

ONE-CARE[®] PRO

PUSH BUTTON SAFETY LANCETS

Ergonomic


Ergonomic T-Shape body for convenient grip

Comfortable

Dual spring safety activation mechanism for enhanced patient's comfort

Easy

Long cap for easy twist off and pull out



PRODUCT ORDERING INFORMATION

Available in six color-coded versions to meet different test blood flow requirements and address a variety of patient skin types.

REF	VERSION	DIAMETER (G) / DEPTH (mm/in)	COLOR	PACKAGING
1011	Micro Flow	30G / 1.2 mm		100 pcs / bx 20 bx (2,000 pcs) / cs
1012	Mini Flow	28G / 1.5 mm		100 pcs / bx 20 bx (2,000 pcs) / cs
1013	Low Flow	25G / 1.6 mm		100 pcs / bx 20 bx (2,000 pcs) / cs
1014	Normal Flow	21G / 1.8 mm		100 pcs / bx 20 bx (2,000 pcs) / cs
1015	High Flow	21G / 2.4 mm		100 pcs / bx 20 bx (2,000 pcs) / cs
1016	Max Flow	21G / 2.8 mm		100 pcs / bx 20 bx (2,000 pcs) / cs



For more information direct
to your smartphone

SCAN THE CODE
> PRO SAFETY LANCETS

ONE-CARE[®]

NEONATAL SAFETY LANCETS



AVAILABLE
1Q-2024

No separate lock or removal of safety tab required, just press the trigger button to activate the device...

One-step and one-hand use for convenience ...

Sterile, single use design with activation indicator for safety ...

Ergonomic shape for convenience ...



Mother with child image for positive feeling, because we care ...

Visual Guides




ORDERING INFORMATION

				
DESCRIPTION	MICRO PREMIE YELLOW	PREMIE PURPLE	NEWBORN GREEN	TODDLER BLUE
INCISION (DEPTH /WIDTH)	0.65 mm / 1.4 mm	0.85 mm / 1.75 mm	1.00 mm / 2.50 mm	2.00 mm / 3.00 mm
INCISION (DEPTH /WIDTH)	REF 1401 50 pcs/box 10 boxes/case	REF 1402 50 pcs/box 10 boxes/case	REF 1403 50 pcs/box 10 boxes/case	REF 1404 50 pcs/box 10 boxes/case

ONE-CARE®

SAFETY BLOOD COLLECTION NEEDLES



Ultra-sharp, lubricated needle for patient comfort

High Fill Rate

Thin-Wall Technology for optimal flow and quick blood collection



Widely accepted safety solution

the side-shield flipped over the needle - minimal training required


Ergonomic design

concave shield to fit thumb shape and textured surface to improve grip

Universal Luer

allows for universal and tight fit to tube holders and multi-sampling

Safety shield activation methods:




- One-hand operation for convenience – thumb or surface activation
- Audible “Click” to indicate that safety shield is activated and the needle is permanently locked



Transparent Tube Holder

- improved visibility of collection tubes filling up with blood,
- tube holders provided separately or pre-attached



ORDERING INFORMATION

CODE	NEEDLE (DIAMETER x LENGTH)	COLOR	PACKAGING
5020	ONE-CARE® Safety Blood Collection Needles - without Holder 22G x 1 1/4"	Black	48 pcs / box 50 bx (2,400 pcs) / case
5021	ONE-CARE® Safety Blood Collection Needles - without Holder 21G x 1 1/4"	Green	48 pcs / box 50 bx (2,400 pcs) / case
5022	ONE-CARE® Safety Blood Collection Needles - without Holder 20G x 1 1/4"	Yellow	48 pcs / box 50 bx (2,400 pcs) / case
5020-H	ONE-CARE® Safety Blood Collection Needles - with Holder 22G x 1 1/4"	Black	100 pcs / box 5 bx (500 pcs) / case
5021-H	ONE-CARE® Safety Blood Collection Needles - with Holder 21G x 1 1/4"	Green	100 pcs / box 5 bx (500 pcs) / case
5022-H	ONE-CARE® Safety Blood Collection Needles - with Holder 20G x 1 1/4"	Yellow	100 pcs / box 5 bx (500 pcs) / case



Blood Collection Safety Needle



Blood Collection Safety Needle with Pre-Attached Holder

ONE-CARE®

SAFETY BLOOD COLLECTION WINGED SETS **GLIDE & LOCK**



Tri-bevel, ultra-sharp, electropolished and lubricated needle for patient comfort

Thin-Wall Technology Canula provides for optimal flow and higher tube fill-rates



Ergonomic wing design permits improved positioning during application

GLIDE & LOCK Safety Shield

enables one-hand activation of safety mechanism, audible click confirms that needle is securely locked



AFTER USE

Transparent Tube Holder

for better visibility of collection tubes filling up with blood, tube holders provided separately or pre-attached

Multi-sampling, universal Luer

allows for tight fit to tube holders and multiple blood sample collection from one patient

DEPF-free, low memory and flexible tubing

assists in handling and preserves the blood chemical and physiological composition

ORDERING INFORMATION

CODE	DESCRIPTION - NEEDLE (DIAMETER x LENGTH)	COLOR	PACKAGING
5040	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Luer-Adapter, 12" Tube	Orange	50 pcs / box 20 bx (1000 pcs) / case
5041	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 24G x 3/4", Luer-Adapter, 12" Tube	Purple	50 pcs / box 20 bx (1000 pcs) / case
5042	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Luer-Adapter, 12" Tube	Blue	50 pcs / box 20 bx (1000 pcs) / case
5043	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 22G x 3/4", Luer-Adapter, 12" Tube	Black	50 pcs / box 20 bx (1000 pcs) / case
5044	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Luer-Adapter, 12" Tube	Green	50 pcs / box 20 bx (1000 pcs) / case
5045	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 20G x 3/4", Luer-Adapter, 12" Tube	Yellow	50 pcs / box 20 bx (1000 pcs) / case
5046	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 19G x 3/4", Luer-Adapter, 12" Tube	Tan	50 pcs / box 20 bx (1000 pcs) / case
5040-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Pre-Attached Holder, 12" Tube	Orange	25 pcs / box 20 bx (500 pcs) / case
5041-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 24G x 3/4", Pre-Attached Holder, 12" Tube	Purple	25 pcs / box 20 bx (500 pcs) / case
5042-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Pre-Attached Holder, 12" Tube	Blue	25 pcs / box 20 bx (500 pcs) / case
5043-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 22G x 3/4", Pre-Attached Holder, 12" Tube	Black	25 pcs / box 20 bx (500 pcs) / case
5044-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Pre-Attached Holder, 12" Tube	Green	25 pcs / box 20 bx (500 pcs) / case
5045-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 20G x 3/4", Pre-Attached Holder, 12" Tube	Yellow	25 pcs / box 20 bx (500 pcs) / case
5046-H	ONE-CARE® GLIDE & LOCK Safety Blood Collection Winged Sets 19G x 3/4", Pre-Attached Holder, 12" Tube	Tan	25 pcs / box 20 bx (500 pcs) / case



ONE-CARE®

SAFETY BLOOD COLLECTION
WINGED SETS **PRESS & LOCK**



Tri-bevel, ultra-sharp, electropolished and lubricated needle
for patient comfort

Thin-Wall Technology Canula provides for optimal flow and higher tube fill-rates



Ergonomic wing design permits improved positioning during application

PRESS & LOCK Safety Shield

enables activation of safety shield in-vein or immediately after the blood draw.

- Press the side buttons and pull back to retract the needle into the shield
- Audible click confirms that needle is securely locked









AFTER USE



Transparent Tube Holder for better visibility of collection tubes filling up with blood, tube holders provided separately or pre-attached

Multi-sampling, universal design Luer allows for tight fit to tube holders and multiple blood sample collection from one patient

DEPF-free, low memory and flexible tubing assists in handling and preserves the blood chemical and physiological composition

-  **1 REMOVE PROTECTIVE CAP**
-  **2 PERFORM VENIPUNCTURE**
-  **3 PRESS SIDE BUTTONS OF THE SAFETY MECHANIS**
-  **4 PULL BACK THE NEEDLE HOLDER AND TUBE UNTIL YOU HEAR "CLICK"**
-  **5 ENSURE THAT NEEDLE IS COMPLETELY LOCKED IN SAFETY SHIELD**
-  **6 DISPOSE OF ALL USED MATERIALS IN A SHARPS CONTAINER**

ORDERING INFORMATION

CODE	DESCRIPTION - NEEDLE (DIAMETER x LENGTH)	COLOR	PACKAGING
5060	ONE-CARE® PRESS & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Luer Adapter, 12" Tube	Orange	50 pcs / box 4 bx (200 pcs) / case
5061	ONE-CARE® PRESS & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Luer Adapter, 12" Tube	Blue	50 pcs / box 4 bx (200 pcs) / case
5062	ONE-CARE® PRESS & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Luer Adapter, 12" Tube	Green	50 pcs / box 4 bx (200 pcs) / case
5060-H	ONE-CARE® PRESS & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Pre-attached Holder, 12" Tube	Orange	25 pcs / box 20 bx (500 pcs) / case
5061-H	ONE-CARE® PRESS & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Pre-attached Holder, 12" Tube	Blue	25 pcs / box 20 bx (500 pcs) / case
5062-H	ONE-CARE® PRESS & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Pre-attached Holder, 12" Tube	Green	25 pcs / box 20 bx (500 pcs) / case



ONE-CARE®

SAFETY WINGED BLOOD COLLECTION SETS

ERGO-SHIELD

● **DEPF-free, low memory and flexible tubing**
assists in handling and preserves the blood chemical and physiological composition

● **Transparent Tube Holder**

for better visibility of collection tubes filling up with blood, tube holders provided separately or pre-attached



● **Tri-bevel, ultra-sharp, electropolished and lubricated needle**
for patient comfort







● **Ergonomic wing design**
permits improved positioning during application



● **Ergonomic Safety Shield**

enables one-hand activation of safety mechanism, audible click confirms that needle is securely locked



Without Holders



REF	PRODUCT	COLOR	PACKAGING
5050	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 25G x 3/4", Luer Adapter, 12" Tube	Orange	100 pcs / box 800 pcs / case
5051	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 23G x 3/4", Luer Adapter, 12" Tube	Blue	100 pcs / box 800 pcs / case
5052	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 21G x 3/4", Luer Adapter, 12" Tube	Green	100 pcs / box 800 pcs / case

With Tube Holders




REF	PRODUCT	COLOR	PACKAGING
5050-H	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 25G x 3/4", Pre-Attached Holder, 12" Tube	Orange	50 pcs / box 400 pcs / case
5051-H	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 23G x 3/4", Pre-Attached Holder, 12" Tube	Blue	50 pcs / box 400 pcs / case
5052-H	ONE-CARE® ERGO-SHIELD Safety Blood Collection Winged Sets 21G x 3/4", Pre-Attached Holder, 12" Tube	Green	50 pcs / box 400 pcs / case

ONE-CARE®

PUSH BUTTON AUTO-RETRACTABLE SAFETY BLOOD COLLECTION WINGED SETS

AVAILABLE
1Q-2024



FEATURES

Provides an ultimate protection of healthcare professionals and patients with the push button and automated needle retraction technology. ONE-CARE PUSH BUTTON Safety Blood Collection Set are sterile, closed systems for venous blood collection. Push-button needle retraction can be activated in-vein with one hand to prevent needle stick injuries.

- **Auto-Retractable Safety Mechanism** - press the button to automatically retract needle into shield
- **One-Hand In-Vein Activation** – in-vein activation with one hand to prevent needle stick injuries
- **Non-curl Tubing** – premium material without “coil memory”
- **Ultra-sharp Needle for Patients’ Comfort**
- **Ergonomic, Enlarged Wings** – secured and convenient handling
- **Enhanced Safety Cap** – improved design of needle safety cap for easy removal
- **High Fill Rate** - Thin-Wall Technology, Lubricated Cannula Needles
- **Regulatory compliant:** Latex & DEHP free materials
- **Cost savings vs major brands:** BD Vacutainer Push Button, Greiner-Bio-One Vacuette Evoprotect

ORDERING INFORMATION

	REF	PRODUCT DESCRIPTION		PACKAGING
<i>Without holders</i>	5070	ONE-CARE® PUSH BUTTON Auto-Retractable Safety Blood Collection Wing Sets 25G x 3/4" Blue, Luer Adapter (without Tube Holder), 12" tube	Blue	50 pcs / bx 4 bx (200 pcs) / cs
	5071	ONE-CARE® PUSH BUTTON Auto-Retractable Safety Blood Collection Wing Sets 23G x 3/4" Blue, Luer Adapter (without Tube Holder), 12" tube	Blue	50 pcs / bx 4 bx (200 pcs) / cs
	5072	ONE-CARE® PUSH BUTTON Auto-Retractable Safety Blood Collection Wing Sets 21G x 3/4" Green, Luer Adapter (without Tube Holder), 12" tube	Green	50 pcs / bx 4 bx (200 pcs) / cs
<i>Pre-attached holders</i>	5070-H	ONE-CARE® PUSH BUTTON Auto-Retractable Safety Blood Collection Wing Sets 25G x 3/4" Blue, Pre-attached Tube Holder, 12" tube	Blue	50 pcs / bx 4 bx (200 pcs) / cs
	5071-H	ONE-CARE® PUSH BUTTON Auto-Retractable Safety Blood Collection Wing Sets 23G x 3/4", Sky Blue, Pre-attached Tube Holder, 12" tube	Blue	50 pcs / bx 4 bx (200 pcs) / cs
	5072-H	ONE-CARE® PUSH BUTTON Auto-Retractable Safety Blood Collection Wing Sets 21G x 3/4", Green, Pre-attached Tube Holder, 12" tube	Green	50 pcs / bx 4 bx (200 pcs) / cs



FEATURES

- Screw-on needle hub compatible with all available holders
- High quality lancet point, silicone treated needle tip for smooth penetration
- Thin wall tube to ensure optimal flow
- Sterile, Non-pyrogenic, non-toxic
- Single use only
- Color coding for easier identification of needle size diameter
- 5 year shelf life, Compliant with ISO 7864:2016

MATERIALS

- ✓ Needle hub - Polypropylene + color masterbatch
- ✓ Cannula – Stainless steel 1.4301 (medical grade)
- ✓ Silicone (medical grade)
- ✓ Glue (medical grade)
- ✓ Rubber sleeve – synthetic rubber (LATEX free)
- ✓ Long protective cap – Polypropylene
- ✓ Short protective cap – Polypropylene



ORDERING INFORMATION

REF	PRODUCT	COLOR	PACKAGING
5001	ONE-CARE® Blood Collection Needles 23G x 1"	Blue	50 pcs / bx 4 bx (200 pcs) / cs
5002	ONE-CARE® Blood Collection Needles 23G x 1 1/2"	Blue	
5003	ONE-CARE® Blood Collection Needles 22G x 1"	Black	
5004	ONE-CARE® Blood Collection Needles 22G x 1 1/2"	Black	
5005	ONE-CARE® Blood Collection Needles 21G x 1"	Green	
5006	ONE-CARE® Blood Collection Needles 21G x 1 1/2"	Green	
5007	ONE-CARE® Blood Collection Needles 20G x 1"	Yellow	
5008	ONE-CARE® Blood Collection Needles 20G x 1 1/2"	Yellow	
5009	ONE-CARE® Blood Collection Needles 18G x 1"	Pink	
5010	ONE-CARE® Blood Collection Needles 18G x 1 1/2"	Pink	



FEATURES

- Size: 0.9x19
- Blood collection needle with holder
- Vacuum Needle
- Instrument Classification Class II

MATERIALS

- ✓ Medical grade PP
- ✓ SUS304
- ✓ ABS

REGULATORY APPROVALS

- US FDA
- CE

STERILIZATION / EXPIRY

- ETO
- Shelf life: 5 years from mfg date

ONE-CARE Blood Collection Winged Sets 25G x 3/4", Luer Adapter without Holder, 12" Tube

ORDERING INFORMATION

	REF	PRODUCT	COLOR	PACKAGING
Blood Collection Wing Sets without Holders	5030	ONE-CARE Blood Collection Winged Sets 25G x 3/4", Luer Adapter without Holder, 12" Tube	Orange	100 pcs / bx 20 bx (2,000 pcs) / cs
	5031	ONE-CARE Blood Collection Winged Sets 24G x 3/4", Luer Adapter without Holder, 12" Tube	Purple	
	5032	ONE-CARE Blood Collection Winged Sets 23G x 3/4", Luer Adapter without Holder, 12" Tube	Light Blue	
	5033	ONE-CARE Blood Collection Winged Sets 22G x 3/4", Luer Adapter without Holder, 12" Tube	Black	
	5034	ONE-CARE Blood Collection Winged Sets 21G x 3/4", Luer Adapter without Holder, 12" Tube	Green	
	5035	ONE-CARE Blood Collection Winged Sets 20G x 3/4", Luer Adapter without Holder, 12" Tube	Yellow	
	5036	ONE-CARE Blood Collection Winged Sets 19G x 3/4", Luer Adapter without Holder, 12" Tube	Brown	
Blood Collection Wing Sets with Holders	5030-H	ONE-CARE Blood Collection Winged Sets 25G x 3/4", Pre-attached Holder, 12" Tube	Orange	100 pcs / bx 20 bx (2,000 pcs) / cs
	5031-H	ONE-CARE Blood Collection Winged Sets 24G x 3/4", Pre-attached Holder, 12" Tube	Purple	
	5032-H	ONE-CARE Blood Collection Winged Sets 23G x 3/4", Pre-attached Holder, 12" Tube	Light Blue	
	5033-H	ONE-CARE Blood Collection Winged Sets 22G x 3/4", Pre-attached Holder, 12" Tube	Black	
	5034-H	ONE-CARE Blood Collection Winged Sets 21G x 3/4", Pre-attached Holder, 12" Tube	Green	
	5035-H	ONE-CARE Blood Collection Winged Sets 20G x 3/4", Pre-attached Holder, 12" Tube	Yellow	
	5036-H	ONE-CARE Blood Collection Winged Sets 19G x 3/4", Pre-attached Holder, 12" Tube	Brown	

ONE-CARE[®]

BLOOD COLLECTION TUBE HOLDERS



Access Holder

Blood Transfer Device



Standard Tube Holder



Safety Holder

Features

- Size: 0.9 x 19 mm
- Transparent material for better visibility of collection tubes filling up with blood
- Universal design for compatibility with blood collection needles and winged sets
- Tight fit connectors to prevent leakage
- DEHP-free & Latex-free materials

Sterilization / Expiry

- ETO
- Shelf life: 5 years from mfg date

Regulatory Compliance

- US FDA, CE, ISO 13485

Materials

- Medical grade PP
- SUS304
- ABS
- All Materials are latex & DEHP free

ORDERING INFORMATION

CODE	DESCRIPTION - NEEDLE (DIAMETER x LENGTH)	COLOR	PACKAGING
5015	ONE-CARE [®] Tube Holder	Transparent	100 pcs / box 20 bx (2000 pcs) / case
5016	ONE-CARE [®] Safety Tube Holder	Transparent/Pink	100 pcs / box 10 bx (1000 pcs) / case
5017	ONE-CARE [®] Blood Transfer Device with Female Luer-Lock and Back-End Needle	Transparent/Pink	200 pcs / box 4 bx (800 pcs) / case
5018	ONE-CARE [®] Access Holder with Male Luer Adapter and Back-End Needle	Transparent/Blue	200 pcs / box 4 bx (800 pcs) / case

WWW.ONE-CARE.COM

ONE-CARE[®] PRO

GLUCOSE MONITOR & STRIPS
FOR PROFESSIONAL USE



Auto
Calibration



0.5 µL Blood
Sample



Alternative
Site
Testing



Instant 14 Days
Average




5 Second
Test



500 Memory
Sets



PC
Download



Strip
Ejector

Bigger Buttons in a compact device

Easy to use



Auto Coding

Simplified way to test,
easy to use, no more
misreading due to
incompatible code card

GOD Enzyme

0.5 µl blood sample
Highly Accurate for
Capillary Whole Blood
samples



For more information direct
to your smartphone

SCAN THE CODE
> GLUCOSE MONITOR &
STRIPS FOR MULTI-PATIENT



Hematocrit | 30-55%
Test Sample | Capillary Whole Blood
Measuring Unit | mg/dL or mmol/L
Measuring Range | 20 - 600 mg/dL (1,1 - 33,3 mmol/L)
Average | 14/07/30 days
Reminder Alarm | 4 sets
External Output | UART with mini USB socket (Optional)
 | Bluetooth (Optional)
Power Supply | 2 AAA batteries
Battery Life | Approximately 1,000 tests
Dimensions | 100 (L) x 57 (L) x 24 (A) mm
Weight | 58 g without batteries

Test Strip Ejector

Reduce the possibility of cross-contamination due to bloodborne pathogens



Test Strip

Longer shelf life 24 months



PRODUCT ORDERING INFORMATION

REF	PRODUCT	PACKAGING	BOXES per CASE
4001	ONE-CARE® PRO GLUCOSE METER FOR PROFESSIONAL USE	1	80
4002	ONE-CARE® PRO GLUCOSE STRIPS FOR PROFESSIONAL USE	50	560
4010	ONE-CARE® CONTROL SOLUTION LEVEL 1 & 3 (LOW & HIGH)	2	900

ONE-CARE®

SAFETY PEN NEEDLES **OPTI**



Safety
Clear safety shield automatically locks after use to protect the needle



Simple
Color indicator becomes invisible after engagement of safety mechanism

Ergonomic Design
Ergonomic main housing with ribs for comfortable grip

Universal
Compatible fit with most pen injectors

UNIVERSAL FIT



Fit all Type-A insulin pen injector as per ISO 11608-2, including:

Compatibility

FlexTouch® Novo Nordisk
FlexPen® Novo Nordisk
NovoPen®Echo Novo Nordisk
NovoPen® 4 Novo Nordisk
NovoPen® 5 Novo Nordisk
Lantus SoloStar®
OptiPen® Pro 1 Sanofi
KwikPen® Humalog

Liprolog® Eli Lilly
Huminsulin® Eli Lilly
HumaPen® Luxura™ Eli Lilly
HumaPen® Luxura™ HD Eli Lilly
HumaPen® MEMOIR™ Eli Lilly
Exenatide Pen Byetta
Autopen® Classic Owen Mumford
Autopen® 24 Owen Mumford

Omnican® Pen 31B.Braun
Ypsopen® Ypsomed
OptiClik® Ypsomed
BerlinPen® Areo Berlin-Chemie
BerlinPen® Areo 2 Berlin Chemie
BerlinPen® 302 Berlin Chemie
Tactipen® Haselmeier
and many others



PRODUCT SELECTION CHART

REF	DIAMETER (G) / DEPTH (mm/in)	COLOR	PACKAGING
3301	30G x 5mm (3/16")	■	100 pcs / bx 20 bx (2,000 pcs) / cs
3302	30G x 8mm (5/16")	■	100 pcs / bx 20 bx (2,000 pcs) / cs

ONE-CARE[®]

SAFETY PEN NEEDLES **PRO**



Insulin delivery has never been easier...

Safety

Clear safety shield automatically locks after use to protect the needle

Ergonomic Design

Ergonomic main housing with ribs for comfortable grip

Simple

Color indicator becomes invisible after engagement of safety mechanism

Universal

Compatible fit with most pen injectors



SAFE

- Automated shielding mechanism ensuring that needle is covered before and after use
- Single-use mechanism preventing sharps injuries and reuse of the device
- All safety pen needles are ETO-sterilized and protected with cover and sterility peel-off tab

SIMPLE

- Safety mechanism engagement indicator via "Click" sound and disappearance of color coded inner shield after fully inserting the needle into the skin
- Optimally designed, clear shield for precise placement on injection site with improved visibility

GENTLE

- Design with concealed needle helps to relieve patient's anxiety
- Polished, silicone-coated, ultra-thin needles are ideally positioned during skin penetration thanks to precise tracks for the safety mechanism elements enhancing patients comfort
- Ribbed housing provides a comfortable grip for easy attachment and removal from the pen

UNIVERSAL

- Universal design for compatible use with most injection pens available on the market
- Tight fit prevents leakages of insulin
- Easy to recognize color-coded product versions address a variety of patient skin types and injection requirements

PRODUCT SELECTION CHART

REF	DESCRIPTION	DIAMETER (G) X DEPTH (MM/IN)	COLOR	PACKAGING
3201	ONE-CARE [®] PRO Safety Pen Needles	31 G x 5mm (3/16")	PURPLE	100 pcs / box 20 bx (2,000 pcs) / cs
3202	ONE-CARE [®] PRO Safety Pen Needles	31 G x 8mm (5/16")	BLUE	100 pcs / box 20 bx (2,000 pcs) / cs

ONE-CARE[®]

SAFETY PEN NEEDLES



Safety

Clear safety shield automatically locks after use to protect the needle

Simple


Color indicator becomes invisible after engagement of safety mechanism

Ergonomic Design

Ergonomic main housing with ribs for comfortable grip

Universal

Compatible fit with most pen injectors



PRODUCT ORDERING INFORMATION

Available in two color-coded product versions for the right fit for every patient and skin type.

REF	VERSION	DIAMETER (G) / DEPTH (mm/in)	COLOR	PACKAGING
3101	Micro	30G x 5mm (3/16")	■	100 pcs / bx 24 bx (2,400 pcs) / cs
3102	Normal	30G x 8mm (5/16")	■	100 pcs / bx 24 bx (2,400 pcs) / cs




For more information direct to your smartphone

SCAN THE CODE
> SAFETY PEN NEEDLES

ONE-CARE[®]

SAFETY PEN NEEDLES **DUO**



Injection needle automatically shielded before & after use










"Thin-Wall Technology" for optimal insulin flow

Ribbed main housing for convenient attachment to pen


Back-end needle automatically covered after removal from pen

FEATURES

- Injection needle automatically shielded before and after use
- Back-end needle automatically shielded after removal of the pen needle from the injection pen
- Glass-clear shield for ease of priming
- Audible "Click" for Activation Indication
- Ribbed main housing for ergonomic design
- Red color visual indication to easily distinguished between used and non-used device
- Tri-bevel, electropolished and siliconized needle for gentle injection
- Thin-Wall technology for optimal insulin flow



ORDERING INFORMATION

REF	NEEDLE	COLOR	PACKAGING
3350	30G x 3/16" 0.3mm x 5mm		100 pcs / bx 8 bx (800 pcs) / cs

ONE-CARE®

SAFETY PEN NEEDLES **PRO DUO**



Safe

Improved visibility during priming, automatic shielding and red color indicator after use

Universal

Compatible fit with most pen injectors

Dual Safety

Back-end needle automatically covered after removal from pen


Ergonomic

Ergonomic main housing for comfortable grip



FEATURES

- Needle automatically shielded before and after use to prevent accidental needle-sticks
- Low resistance during injection and activation of the safety mechanism
- Red-color activation indicator to differentiate before used and non-used device
- See-through window for improved visibility during priming
- Ergonomic design for easy attachment into pens
- Optimal skin contact platform that promotes 90-degree injection angle
- Compatible with most popular pens on the market
- Different needle configurations available



ORDERING INFORMATION

REF	NEEDLE	COLOR	PACKAGING
3360	ONE-CARE PRO DUO Safety Pen Needles 31G x 5mm	■	100 pcs / box 1200 pcs / case
3361	ONE-CARE PRO DUO Safety Pen Needles 31G x 8mm	■	100 pcs / box 1200 pcs / case

ONE-CARE[®]

INSULIN AND TB SYRINGES

Gentle

Precisely ground, polished, and silicone-coated needles for patient's comfort

Safety Shield

Clear safety shield easy to push forward for transportation or push and twist to permanently lock for safety

Clear, easy to read marking

Clear marking and smooth movement of plunger for precise dosing.

No medication loss

Minimal dead space to avoid medication loss



PRODUCT ORDERING INFORMATION

	CODE	SIZE	NEEDLE (DIAMETER x LENGTH)	PACKAGING
INSULIN SAFETY SYRINGES	2101	1 mL	29G x 1/2" 0.3mm x 12.7mm	100 pcs / bx 10 bx (1,000 pcs) / cs
	2102	0.5 mL	29G x 1/2" 0.3mm x 12.7mm	100 pcs / bx 10 bx (1,000 pcs) / cs

	CODE	SIZE	NEEDLE (DIAMETER x LENGTH)	PACKAGING
TUBERCULIN SAFETY SYRINGES	2107	1 mL	28G x 1/2" 0.36mm x 12.7mm	100 pcs / bx 10 bx (1,000 pcs) / cs
	2108	1 mL	27G x 1/2" 0.40mm x 12.7mm	100 pcs / bx 10 bx (1,000 pcs) / cs
	2109	1 mL	25G x 5/8" 0.50mm x 16.0mm	100 pcs / bx 10 bx (1,000 pcs) / cs



SCAN THE CODE
> INSULIN AND TUBERCULIN SAFETY SYRINGES

For more information direct to your smartphone

ONE-CARE[®]









INSULIN SAFETY SYRINGES **OPTI**



- Innovative Safety Shield
- One-Step Activation
- One-Hand Operation
- Low Dead Space



Push forward the shield with thumb to cover the needles and lock the device



ORDERING INFORMATION


REF	SYRINGE SIZE	NEEDLE	COLOR	PACKAGING
2141	1 mL	29G x 1/2" 0.33 x 12.7 mm	■	100 pcs / bx 8 bx (800 pcs) / cs
2142	0.5 mL	29G x 1/2" 0.33 x 12.7 mm	■	100 pcs / bx 8 bx (800 pcs) / cs

ONE-CARE®

AUTO-RETRACTABLE
INSULIN SAFETY SYRINGES




**AVAILABLE
1Q-2024**



FEATURES

Provides an ultimate protection of healthcare professionals and patients with automated needle retraction technology. Safety mechanism of ONE-CARE Auto-Retractable safety syringes can be activated with one hand while the needle is in skin to prevent needle stick injuries.

- **Market proven** technique of use
- **One-hand operation**
- **Patented** claw design to prevent falling before using
- **Precise needle point** design and grinding, silicone-coated for smooth insertion
- **Latex and DEHP free materials** for regulatory compliance
- **Clear graduation** for accuracy in small dosage injection
- **27-31G sizes available**



ORDERING INFORMATION

REF	PRODUCT DESCRIPTION	PACKAGING
2211	ONE-CARE® Auto-Retractable Insulin Safety Syringe 1ml 31G x 5/16"	100 pcs / bx 12 bx (1,200 pcs) / cs
2212	ONE-CARE® Auto-Retractable Insulin Safety Syringe 1ml 31G x 1/2"	100 pcs / bx 12 bx (1,200 pcs) / cs
2213	ONE-CARE® Auto-Retractable Insulin Safety Syringe 1ml 30G x 1/2"	100 pcs / bx 12 bx (1,200 pcs) / cs
2214	ONE-CARE® Auto-Retractable Insulin Safety Syringe 1ml 29G x 1/2"	100 pcs / bx 12 bx (1,200 pcs) / cs
2215	ONE-CARE® Auto-Retractable Insulin Safety Syringe 0.5ml 29G x 1/2"	100 pcs / bx 12 bx (1,200 pcs) / cs
2216	ONE-CARE® Auto-Retractable Insulin Safety Syringe 1ml 27G x 1/2"	100 pcs / bx 12 bx (1,200 pcs) / cs

Safety Syringes and Safety Needles

Luer-Lock Syringes 1ml with Safety Needles 23G x 1" (25mm)



REF 2110

1 ML 23G x 1" (25mm)

Easy, one-hand activation of safety shield

Audible click indicates safety mechanism activation

Locking mechanism ensures needle is secure and completely locked in the safety shield

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	PRODUCT / APPLICATION	PACKAGING
2110	ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 23G x 1"(25mm)	Syringes with Safety Needles for Injections	100 pcs / bx 12 bx / cs

Luer-Lock Syringes 1ml with Safety Needles 25G x 1" (25mm)



REF 2112

1 ML 25G x 1" (25mm)

Easy, one-hand activation of safety shield

Audible click indicates safety mechanism activation

Locking mechanism ensures needle is secure and completely locked in the safety shield

ORDERING INFORMATION

REF	PRODUCT DESCRIPTION	PRODUCT / APPLICATION	PACKAGING
2112	ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 25G x 1"(25mm)	Syringes with Safety Needles for Injections	100 pcs / bx 12 bx / cs

Safety Syringes and Safety Needles

SAFETY NEEDLES

ONE-CARE® Safety Needles

Safety Needles 23Gx1" (25mm)

REF 2111



Safety Needles 25Gx1" (25mm)

REF 2113



Safety Needles 21Gx1.5" (38mm)

REF 2114




Safety Needles 18Gx1.5" (38mm)


REF 2115



Before Use



After Use



ORDERING INFORMATION

REF	PRODUCT DESCRIPTION	PRODUCT / APPLICATION	PACKAGING
2111	ONE-CARE® Safety Needles 23G x 1"(25mm)	Safety Needles	100 pcs / bx 20 bx / cs
2113	ONE-CARE® Safety Needles 25G x 1"(25mm)	Safety Needles	100 pcs / bx 20 bx / cs
2114	ONE-CARE® Safety Needles 21G x 1.5"(38mm)	Safety Needles	100 pcs / bx 20 bx / cs
2115	ONE-CARE® Safety Needles 18G x 1.5"(38mm)	Safety Needles	100 pcs / bx 20 bx / cs
2116	ONE-CARE® Safety Needles 27G x 3/4"(16mm)	Safety Needles	100 pcs / bx 20 bx / cs

ONE-CARE®

Syringes and Needles

LUER-SLIP DISPOSABLE SYRINGES, 3-PART

REF 2001

1 ml



Specially designed 1 ml syringes for accurate dose administration

Silicone lubrication

Ensures smooth plunger movement

Safe plunger backstop

Prevents accidental withdrawal of plunger

Large and stable plunger grip

Ensures stability and comfort during injection

High-transparency barrel

Makes it easy to see syringe content

PRODUCT/APPLICATION	PRODUCT DESCRIPTION	REF	BOX / CASE
Disposable Syringes for Vaccine Injections	ONE-CARE® Luer-Slip Syringes 1ml, 3-part	2001	100 pcs / bx 24 bx / cs

LUER-LOCK DISPOSABLE SYRINGES, 3-PART

REF 2002 - 2003 - 2004 - 2005

1 mL

3 mL

5 mL

10 mL



Luer-Lock 1ml syringes for secure connection with hypodermic needles

Silicone lubrication

Ensures smooth plunger movement

High-transparency barrel

Makes it easy to see syringe content

Safe plunger backstop

Prevents accidental withdrawal of plunger

Large and stable plunger grip

Ensures stability and comfort during injection

PRODUCT/APPLICATION	PRODUCT DESCRIPTION	REF	BOX / CASE
Disposable Syringes for Vaccine Injections	ONE-CARE® Luer-Lock Syringes 1mL, 3-part	2002	100 pcs / box 24 bx / cs
	ONE-CARE® Luer-Lock Syringes 3mL, 3-part	2003	100 pcs / box 32 bx / cs
	ONE-CARE® Luer-Lock Syringes 5mL, 3-part	2004	100 pcs / box 16 bx / cs
	ONE-CARE® Luer-Lock Syringes 10mL, 3-part	2005	100 pcs / box 15 bx / cs

ONE-CARE®

Syringes and Needles

DISPOSABLE NEEDLES

ONE-CARE® DISPOSABLE NEEDLES 23G X 1" (25mm) REF 23100-CM

ONE-CARE® DISPOSABLE NEEDLES 23G X 1.5" (38mm) REF 23150-CM



ONE-CARE® hypodermic needles for vaccine administration have precise tri-bevel design, are electro-polished and lubricated for gentle and virtually pain-free injections.

ONE-CARE® DISPOSABLE NEEDLES 21G X 1.5" (38mm) REF 21150-CM










ONE-CARE® Canula Needles 21G x 1.5" (38mm) are recommended for vaccine preparation as may be required by vaccine manufacturer's instructions for use.

PRODUCT/APPLICATION	PRODUCT DESCRIPTION	REF	BOX / CASE
Hypodermic needles for vaccine injections	ONE-CARE® Hypodermic Needles 23G x 1.5" (38mm)	23150-CM	100 pcs / bx 100 bx (10,000 pcs) / cs
	ONE-CARE® Hypodermic Needles 23G x 1" (25mm)	23100-CM	
	ONE-CARE® Hypodermic Needles 23G x 5/8" (16mm)	23058-CM	
	ONE-CARE® Hypodermic Needles 25G x 1.5" (38mm)	25150-CM	
	ONE-CARE® Hypodermic Needles 25G x 1" (25mm)	25100-CM	
	ONE-CARE® Hypodermic Needles 25G x 5/8" (16mm)	25058-CM	

ONE-CARE®

GLUCOSE MONITORING SYSTEM
FOR SELF-TESTING



- 
FREE
 Auto Calibration
- 
0.5µl
 0.5 µL Blood Sample
- AST**
 Alternative Site Testing
- 
5^{sec}
 5 Second Test
- 
14
 Instant 14 Days Average
- 
500
 500 Memory Sets
- 
PC Download
- 
Strip Ejector

TEST KIT

- > 1 Blood Glucose Meter with batteries
- > 10 Blood Glucose Test Strips
- > 1 Adjustable Depth Lancing Device
- > 10 Sterilized Lancets
- > 1 Protective Carrying Case
- > User Manual, Package Inserts

Note: The contents of a Complete Kit may vary and are subject to change according to the actual purchase.



PRODUCT ORDERING INFORMATION

REF	VERSION	PCS per BOX	BOXES per CASE
4006	ONE-CARE® GLUCOSE MONITORING KIT FOR SELF-TESTING	1	80
4007	ONE-CARE® GLUCOSE STRIPS FOR SELF-TESTING	50	560



For more information direct to your smartphone

SCAN THE CODE
> GLUCOSE METER

ONE-CARE[®]

PEN NEEDLES



Gentle

Precisely grinded three bevel tip, electro-polished and lubricated canulas for enhanced comfort

Optimal Insulin Flow

Thin- and Extra-Thin-Wall Technology

Universal Fit

ISO 11608-2 compliant and universal thread design for compatibility with all type A insulin pen injectors

PRODUCT ORDERING INFORMATION

Available in four color-coded versions to meet different test blood flow requirements and address a variety of patient skin types.



REF	DIAMETER (G) x DEPTH (mm/in)	COLOR	PACKAGING
3011	33G x 4mm	Pink	100 pcs / bx 50 bx (5,000 pcs) / cs
3012	32G x 4mm	Light Green	100 pcs / bx 50 bx (5,000 pcs) / cs
3013	32G x 5mm	Brown	100 pcs / bx 50 bx (5,000 pcs) / cs
3014	32G x 6mm	Yellow	100 pcs / bx 50 bx (5,000 pcs) / cs
3015	32G x 8mm	Grey	100 pcs / bx 50 bx (5,000 pcs) / cs
3016	31G x 4mm	Light Blue	100 pcs / bx 50 bx (5,000 pcs) / cs
3017	31G x 5mm	Purple	100 pcs / bx 50 bx (5,000 pcs) / cs
3018	31G x 6mm	Dark Blue	100 pcs / bx 50 bx (5,000 pcs) / cs
3019	31G x 8mm	Orange	100 pcs / bx 50 bx (5,000 pcs) / cs



For more information direct to your smartphone

SCAN THE CODE
> PEN NEEDLES

ONE-CARE[®]

Lancing Device & Lancets



Easy
Easy Depth Setting
with clear marking


"Ready to Activate"
Visual Indicator

Compact Size
Compact Size for discrete
blood sampling

Convenient
Convenient Loading System

Universal
Universal Design for
fit with most common lancets

**Ultra-thin, electro-polished
lancets for comfortable sampling**



**Automated Ejector
for safe lancet removal**



PRODUCT ORDERING INFORMATION

REF	DIAMETER (G) / DEPTH (mm/in)	DESCRIPTION	COLOR	PACKAGING
1200	10 - Depth Setting	Universal Lancing Device	□	100 pcs / bx 100 bx (10,000 pcs) / cs
1201	23G	23G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs
1202	26G	26G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs
1203	28G	28G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs
1204	30G	30G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs
1205	31G	23G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs
1206	32G	32G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs
1207	33G	33G Universal Blood Lancets	■	100 pcs / bx 100 bx (10,000 pcs) / cs



For more information direct
to your smartphone

SCAN THE CODE
> LANCING DEVICE & LANCETS

ONE-CARE[®]

ALCOHOL PREP PADS



ONE-CARE™ Alcohol Prep Pads are saturated with 70% Isopropyl Alcohol for optimal antiseptic action (anti-bacteria, anti-virus). They are intended for topical use in preparation of skin prior to an injection or blood collection, cleaning of minor cuts and scrapes, or cleaning of medical equipment and supplies. Reliable 4-layer packaging provides air-tight seal and prevents leakage or dry-out.

FEATURES

- Soft, absorbent, nonwoven pads saturated with 70% v/v Isopropyl Alcohol
- For Individual, Professional and Hospital Use
- Reliable 4-layer packaging
- Compact Box of 100 pads
- Attractive Price

BENEFITS

- Exceeds CDC's recommended 60% alcohol content. Optimal antiseptic action helps to reduce bacteria.
- Intended for use in preparation of skin prior to an injection or blood collection, cleaning of minor cuts and scrapes, or cleaning of medical equipment and supplies
- Air-tight seal prevents leakage or dry-out
- Safe to use in medical setting
- Lower cost of warehousing and shipping
Financial value / cost savings for quality and reliable product as compared to major national brands

PRODUCT ORDERING INFORMATION

REF	DESCRIPTION 1	DESCRIPTION 2	SIZE	PACKAGING
1020	ONE-CARE™ Alcohol Prep Pads	Saturated with 70% Isopropyl Alcohol	M	100 pcs / bx 50 bx (5,000 pcs) / cs
1022	ONE-CARE™ Alcohol Prep Pads	Saturated with 70% Isopropyl Alcohol	M	200 pcs / bx 50 bx (10,000 pcs) / cs



For more information direct to your smartphone

SCAN THE CODE
> ALCOHOL PREP PADS

Nitrile Examination Gloves



MEDICAL GRADE, POWDER-FREE, TEXTURED EXAMINATION GLOVES



Product Description

Type of the glove

Non-sterile, powder free, examination and protective glove for single use

Material

Nitrile

Colour

Blue

Shape

Ambidextrous

Cuff

Beaded

External surface

Finger textured 35 mm. From Tip

Internal surface

Polymerized + chlorinated

Packaging

10 x 100 pcs

Physical Properties

DIMENSIONS

WIDTH (in)
(mm)
Tolerances

S	M	L	XL
3.14 ± 0.39 80 ± 10 (6-7)	3.74 ± 0.39 95 ± 10 (7-8)	4.33 ± 0.39 110 ± 10 (8-9)	>4.33 ± 0.39 >110 ± 10 (9-10)

LENGTH (in)
(mm)
minimum

≥ 9.44 ≥ 240	≥ 9.44 ≥ 240	≥ 9.44 ≥ 240	≥ 9.44 ≥ 240
-----------------	-----------------	-----------------	-----------------

THICKNESS (mm)/(mil)
(Single wall)
minimum

M - Finger: min 3.15 mil | min 0.08 mm
M - Palm: min 2.36 mil | min 0.06 mm
M - Cuff: min 1.97 mil | min 0.05 mm

ELONGATION AT BREAK (%)
minimum

M - 500 (ISO 11193-1)
M - 400 (ISO 11193-1)

FORCE AT BREAK (N)
minimum

6.0 (EN 455-2) / 7.0 (ISO 11193-1)

Manufacturing and Safety standards


AQL	Manufacturing final release: G-I inspection level AQL 1.5 in accordance with ISO 2859-1	
US FDA Registration	Class 1, Product Code LZA, 510(k) K172930, device listing D320597	
CE classification	Class I - Medical Device Regulation (EU) 2017/745	Category III – Personal Protective Equipment (Regulation (EU)
Compliances	EN 455-1, EN 455-2, EN 455-3, EN 455-4 EN ISO 15223-1 EN 1041 ISO 11193-1 EN ISO 13485	EN ISO 374-1 (Type B), EN 374-2, EN 374-4, EN ISO 374-5 EN 16523-1 EN 420
Viral test	Test in accordance with ASTM F1671 & ISO 16604	
Cytostatics permeation	Test in accordance with ASTM D6978	
Chemical substances permeation	Test in accordance with EN 16523-1	
Food contact	Declaration of conformity for food contact in accordance with Regulation (EC) No 1935/2004 and Overall Migration Test in accordance with European Commission Directive 93/11/EEC and Council of Europe Resolution AP (2004) 4	
Shelf life	3 years	
Storage Instruction	Keep out of direct sunlight. Store in a cool, dry place in temperature 5-35° C. Keep away from sources of ozone and ignition	

PRODUCT/APPLICATION	PRODUCT DESCRIPTION	REF	BOX / CASE
Personal Protection	NITRYLEX Nitrile Examination Gloves, Powder-Free, Blue	S: RD30019902 M: RD30019903 L: RD30019904 XL: RD30019905	100 pcs / bx 10 bx (1,000 pcs) / cs

For more details, contact our sales department

Particulate Respirators Foldable-Design Type N95

N95 FOLDABLE-DESIGN RESPIRATOR, NIOSH CERTIFIED, CDC and US FDA Listed



- Soft smooth lining
- Considerable fit design
- adjustable nosepiece
- Comfortable nose foam
- Easy breathing speaking

MS8225 Foldable-Design N95

Mask Type	Standards	Filtration Efficiency *
Respirator Mask (PFE Standard)	USA: NIOSH (42 CFR 84)	N95 PFE ≥ 95% for 0.3-micron particles

* (PFE - Particle Filtration or BFE - Bacteria Filtration Efficiency)

Respiratory Mask of Class N95 is suitable for filtering of the most penetrating particles of 0.3-micron size that can be found in the air.

- 0.1 micron:** Particle Filtration Efficiency standard (PFE).
- 0.3 micron:** Used to represent the most-penetrating particle size (MPPS), which is the most difficult size particle to capture.
- 3.0 microns:** Bacteria Filtration Efficiency standard (BFE) for masks per EN 14683.


Name	MS8225 Folded N95
Colour	White
Item Number	8225
Weight/unit (oz) (g)	±0.34oz 9.6g
Outer layer	PP needle-punched nonwoven
Marking	As per NIOSH
Soft foam nosepiece	Black rectangular foam
Dimensions (in)	4.53*6.30 in (± 0.20in)
Dimensions (mm)	115*160 mm (± 5mm)
Description	N95 Particulate Respirator
Material	4-plys, PP spunbond outer layer, PP meltbown high filtration layer, PP needle-punched nonwoven inner layer
Elastic head band	Spandex headstraps
Inner layer	PP needle-punched nonwoven
Filter media	PP meltbown high filtration layer
Adjustable nose clip	Iron wire + plastic nose clip
Characteristics	Conforms to NIOSH 42 CFR standard. TC#84A-7447 95% filtration effective against particulate aerosols free of oils. Adjustable plastic nose clip, Latex-Free elastic ultrasonic welded to cup

ORDERING INFORMATION	Materials	Boxes	Cases
		350 gsm white cardboard	Brown 5-ply single wall
2225	Dimensions (in) (mm)	5.12*4.72*6.30 (in) 130*120*160 (mm)	24.41*11.02*13.38 (in) 620*280*340 (mm)
	Gross Weight: (lb) 11.90lb/cs (kg) 5.4kg/cs		

Package 20 masks/polybag/box * 20 boxes = 400 masks / case

Particulate Respirators Cone-Design Type N95


N95 CONE-DESIGN DISPOSABLE RESPIRATOR, NIOSH CERTIFIED, CDC and US FDA Listed



- Contour Fit Design
- Adjustable Nosepiece
- Comfortable Nose Foam
- Easy Breathing Speaking
- NIOSH N95 APPROVED

MS6115L Cone-Design Type N95

Name	MS6115L Cup-type N95
Colour	White
Item Number	6115L
Weight/unit (oz) (g)	±0.3oz 8.68g
Outer layer	PP Spunbond
Marking	As per NIOSH
Soft foam nosepiece	Black PU foam
Dimensions (in)	4.72*5.51*2.09 in (±0.20in)
Dimensions (mm)	120*140*53 mm (± 5mm)
Description	N95 Particulate Respirator
Material	3-plys, PP spunbond outer layer, PP meltbown high filtration layer, PP needle-punched nonwoven inner layer
Elastic head band	Latex-free elastic-width 7mm
Inner layer	PP needle-punched nonwoven
Filter media	PP meltbown high filtration layer
Adjustable nose clip	Aluminium
Characteristics	2 CFR standard. TC#84A-5530 95% filtration effective against particulate aerosols, free of oils, Aluminum nose-clip, Latex-Free elastic ultrasonic welded to cup

Mask Type	Standards	Filtration Efficiency*
Respirator Mask (PFE Standard)		N95
	USA: NIOSH (42 CFR 84)	PFE ≥ 95% for 0.3-micron particles

* (PFE - Particle Filtration or BFE - Bacteria Filtration Efficiency)

Respiratory Mask of Class N95 is suitable for filtering of the most penetrating particles of 0.3-micron size that can be found in the air.

- 0.1 micron:** Particle Filtration Efficiency standard (PFE).
- 0.3 micron:** Used to represent the most-penetrating particle size (MPPS), which is the most difficult size particle to capture.
- 3.0 microns:** Bacteria Filtration Efficiency standard (BFE) for masks per EN 14683.

ORDERING INFORMATION	Materials	Boxes	Cases
		350 gsm white cardboard	Brown 5-ply single wall
6115L	Dimensions (in) (mm)	5.31*4.92*7.87 (in) 135*125*200 (mm)	25.59*11.42*16.53 (in) 650*290*420 (mm)
		Gross Weight: (lb) 12.46lb/cs	(kg) 5.65kg/cs

Package 20 masks/polybag/box * 20 boxes = 400 masks / case

ONE-CARE[®]

COLLECT BLOOD SAMPLE

Blood Collection Needles & Winged Sets

Safety Lancets

Consumer Lancets & Lancing Device

PERFORM TEST

Glucose Meters

Rapid Tests

INJECT MEDICATION/VACCINATION

Pen Needles

Insulin & TB Syringes

Syringes & Safety Needles

PROTECT HEALTHCARE PROFESSIONALS AND PATIENTS

Alcohol Prep Pads

Nitrile Exam Gloves

N95 Surgical Masks

www.one-care.com

MediVena[®]

inspired by innovation

AMERICAS:

MediVena

1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA

+1 404 514 2586 info@medivena.com

CUSTOMER SERVICE TOLL FREE: 888 270 6457

EUROPE:

MediVena Sp. z o.o.

Al. Jana Pawla II 27, 00-867 Warsaw, POLAND

+48 736 667 468

www.medivena.com