

ONE-CARE[®]

GLUCOSE MONITORING SYSTEM
FOR SELF-TESTING

Glucose monitoring has never been easier...



 Auto Calibration	 0.5 µL Blood Sample	 Alternative Site Testing	 Instant 14 Days Average
 5 Second Test	 500 Memory Sets	 PC Download	 Strip Ejector



Bigger Buttons in a compact device

Easy to use for elderly people

USB and Bluetooth (Optional)

Provide both wire and wireless connection for data transmission

Test Strip Ejector

Reduce the possibility of cross-contamination due to bloodborne pathogens

Test Strip

Longer shelf life 24 months

Certification



SPECIFICATIONS




Auto Coding

Simplified way to test easy to use, no more misreading due to incompatible code card

GDO Enzyme


0,5 µl blood sample
Accurate for Fresh Capillary Whole Blood samples



Hematocrit	30-55%
Test Sample	Capillary Whole Blood
Measuring Unit	mg/dL or mmol/L
Measuring Range	20 - 600 mg/dL (1,1 - 33,3 mmol/L)
Average	14/07/30 days
Reminder Alarm	4 sets
External Output	UART with mini USB socket (Optional) Bluetooth (Optional)
Power Supply	2 AAA batteries
Battery Life	Approximately 1,000 tests
Dimensions	100 (L) x 57 (L) x 24 (A) mm
Weight	58 g without batteries


INSTRUCTIONS OF USE

1 STEP




Insert a Test Strip

2 STEP




Apply Blood Sample

3 STEP



Wait 5 Seconds and Check Results

4 STEP



Discard the used Test Strip

ORDERING INFORMATION

TEST KIT



TEST KIT

- > 1 Blood Glucose Meter with batteries
- > 10 Blood Glucose Test Strips
- > 1 Adjustable Depth Lancing Device
- > 10 Sterilized Lancets
- > 1 Protective Carrying Case
- > User Manual, Package Inserts

Note: The contents of a Complete Kit may vary and are subject to change according to the actual purchase.

CODE	PRODUCT	PCS per BOX	BOXES per CASE
4006	ONE-CARE® PRO GLUCOSE MONITORING KIT FOR SELF-TESTING	1	80
4007	ONE-CARE® PRO GLUCOSE STRIPS FOR SELF-TESTING	50	560

WWW.ONE-CARE.COM

MediVena
inspired by innovation

Americas

USA
MediVena
1 Glenlake Parkway, Ste 700
Atlanta, GA 30328, USA
+1 404 422 1513
info@medivena.com
www.medivena.com

Europe

Poland
MediVena Sp. z o.o.
Al. Jana Pawla II 27
00-867 Warsaw, Poland
+48 22 299 44 19
info@medivena.com
www.medivena.com

