

ONE-CARE™ PEN NEEDLES



Insulin delivery has never been easier...



Gentle

Precisely grinded three bevel tip, electro-polished and lubricated canulas for enhanced comfort



Universal Fit

for compatibility with most pen injectors

Optimal Insulin Flow

Thin and Extra-Thin-wall technology for fast insulin delivery

GENTLE

- Electro-polished, silicone-coated, ultra-thin needles guarantee enhanced user comfort
- Easy flow cannula minimizes time of insulin delivery
- Smooth and quick insertion thanks to special lubrication makes the injection a barely-noticeable experience

CONVENIENT

- Ergonomic design with ribbed housing ensures comfortable grip for secure attachment onto the pen injector
- Color-coded versions the right needle selection for every patient
- ETO-sterilized and protected with outer cover and easy peel-off tab

UNIVERSAL

- **Twist-on** universal design for compatibility with most injection pens available on the market
- Tight fit prevents leakages and medication loss
- 100% automated camera-control for quality
- Rigorously tested per international standards for fit, functionality and reliable performance

FlexTouch® Novo Nordisk
 FlexPen® Novo Nordisk
 NovoPen® Echo Novo Nordisk
 NovoPen® 4 Novo Nordisk
 NovoPen® 5 Novo Nordisk
 Lantus SoloStar®
 OptiPen® Pro 1 Sanofi
 KwikPen® Humalog


Liprolog® Eli Lilly
 Huminsulin® Eli Lilly
 HumaPen® Luxura™ Eli Lilly
 HumaPen® Luxura™ HD Eli Lilly
 HumaPen® MEMOIR™ Eli Lilly
 Exenatide Pen Byetta
 Autopen® Classic Owen Mumford
 Autopen® 24 Owen Mumford

Omnican® Pen 31B. Braun
 Ypsopen® Ypsomed
 OptiClik® Ypsomed
 BerlinPen® Areo Berlin-Chemie
 BerlinPen® Areo 2 Berlin Chemie
 BerlinPen® 302 Berlin Chemie
 Tactipen® Haselmeier
 and many others




INSTRUCTIONS FOR USE

Intuitive and simple handling. Insulin delivery has never been easier...



Protective brand label



Outer needle cap



Inner needle cap



Hub needle

PRODUCT SELECTION CHART

Available in six color-coded versions to meet different test blood flow requirements and address a variety of patient skin types.

CODE	VERSION	DIAMETER (G) x DEPTH (mm/in)	PACKAGING
3301	MICRO	32G x 4mm (3/16in)	100 pcs / box 2000 pcs / case
3302	MINI	31G x 5mm (5/16in)	100 pcs / box 2000 pcs / case
3303	NORMAL	31G x 6mm (5/16in)	100 pcs / box 2000 pcs / case
3304	MAX	31G x 8mm (5/16in)	100 pcs / box 2000 pcs / case

WWW.ONE-CARE.COM

MediVena
 inspired by innovation

Americas

USA
 MediVena
 1 Glenlake Parkway, Ste 700
 Atlanta, GA 30328, USA
 +1 404 422 1513
 info@medivena.com
 www.medivena.com

Europe

Poland
 MediVena Sp. z o.o.
 Al. Jana Pawla II 27
 00-867 Warsaw, Poland
 +48 22 299 44 19
 info@medivena.com
 www.medivena.com

