# **ONE-CARE**® SYRINGES AND NEEDLES FOR VACCINATION

In response to COVID-19, MediVena is offering **ONE-CARE® Syringes and Needles** that address requirements for sample preparation and vaccine injections. In addition, **ONE-CARE® Alcohol Prep Pads** are provided for disinfection of the surface of the vials and the injection sites.

**ONE-CARE**<sup>®</sup> product configurations recommended in this brochure are based on the vaccination guidelines and instructions for use of the first COVID-19 vaccines to cover all needs of the vaccination programs in different countries.



MediVena's **ONE-CARE®** Syringes and Needles have "low-dead space" designs to reduce waste of vaccination and to allow for extractions as many doses as possible to maximize the vaccination opportunities for population. In particular, **ONE-CARE®** injection syringes of 1 ml have the elongated piston tip design to leave low-to-zero "dead space" within syringe.



## SYRINGES AND NEEDLES FOR VACCINATION

### **ONE-CARE<sup>®</sup> LUER-SLIP DISPOSABLE SYRINGES**

REF 2001



## Specially designed 1 ml syringes for accurate dose administration

**Silicone lubrication** Ensures smooth plunger movement

**Safe plunger backstop** Prevents accidental withdrawal of plunger

Large and stable plunger grip Ensures stability and comfort during injection

**High-transparency barrel** Makes it easy to see syringe content

Luer-Lock 1ml syringes for secure connection with hypodermic needles

Silicone lubrication Ensures smooth plunger movement

High-transparency barrel Makes it easy to see syringe content

Safe plunger backstop Prevents accidental withdrawal of plunger

Large and stable plunger grip Ensures stability and comfort during injection

#### ONE-CARE<sup>®</sup> LUER-LOCK DISPOSABLE SYRINGES REF 2002



**ONE-CARE® LUER-SLIP DISPOSABLE SYRINGES** 

Bill Chart and a state and a state of the st

2ml Luer-Slip 2-part syringe for easy and accurate vaccine preparation

Clear and precise marking for accurate dosing

Sterile, blister-packed

SYRINGES AND NEEDLES FOR VACCINATION



REF 2003

**2 ml** 

## ONE-CARE

## **PRODUCT DESCRIPTION**

## **ONE-**CARE

## ONE-CARE<sup>®</sup> DISPOSABLE NEEDLES 23G X 1" (25mm) REF 23100-CM ONE-CARE<sup>®</sup> DISPOSABLE NEEDLES 23G X 1.5" (38mm) REF 23150-CM



**ONE-CARE®** hypodermic needles for vaccine administration have precise tri-bevel design, are electro-polished and lubricated for gentle and virtually pain-free injections.

#### ONE-CARE® DISPOSABLE NEEDLES 21G X 1.5" (38mm) REF 21150-CM



**ONE-CARE® Canula Needles 21G x 1.5" (38mm)** are recommended for vacinne preparation as may be required by vaccine manufacturer's instructions for use.

#### **ONE-CARE® ALCOHOL PREP PADS**

REF 1021



70% Isopropyl Alcohol For External Antiseptic Use Only Soft, absorbent, nonwoven pads are saturated with 70% Isopropyl Alcohol which exceeds CDC's recommended 60% alcohol content for optimal antiseptic action.

Intended for Individual, Professional and Hospital Use in preparation of skin prior to an injection or blood collection, cleaning of minor cuts and scrapes, or cleaning of medical equipment and supplies

Compact Box of 100 pads for optimized storage and shipping cost

SYRINGES AND NEEDLES FOR VACCINATION



## **ONE-**CARE<sup>®</sup>

| PRODUCT/APPLICATION                              | PRODUCT DESCRIPTION  | REF      | BOX / CASE                  |
|--|--|----------|-----------------------------|
| Disposable<br>Syringes for<br>Vaccine Injections | ONE-CARE <sup>®</sup> Luer-Slip Syringes 1ml, 3-part       | 2001     | 100 pcs / bx<br>24 bx / cs  |
|  | ONE-CARE <sup>®</sup> Luer-Lock Syringes 1ml, 3-part       | 2002     |                             |
| Hypodermic Needles<br>for Vaccine Injections     | ONE-CARE® Hypodermic Needles 23G x 1.5" (38mm)             | 23150-CM | 100 pcs / bx<br>100 bx / cs |
|  | ONE-CARE <sup>®</sup> Hypodermic Needles 23G x 1"(25mm)    | 23100-CM |                             |
|  | ONE-CARE® Hypodermic Needles 23G x 5/8"(16mm)              | 23058-CM |                             |
|  | ONE-CARE <sup>®</sup> Hypodermic Needles 25G x 1.5" (38mm) | 25150-CM |                             |
|  | ONE-CARE <sup>®</sup> Hypodermic Needles 25G x 1.5" (38mm) | 25100-CM |                             |
|  | ONE-CARE <sup>®</sup> Hypodermic Needles 25G x 5/8"(16mm)  | 25058-CM |                             |
| Disposable Syringes<br>for Vaccine Dilution      | ONE-CARE <sup>®</sup> Luer-Slip Syringes 2ml, 2-part       | 2003     | 100 pcs / bx<br>25 bx / cs  |
|  | ONE-CARE <sup>®</sup> Luer-Slip Syringes 3ml, 2-part       | 2004     | 100 pcs / bx<br>25 bx / cs  |
|  | ONE-CARE <sup>®</sup> Luer-Slip Syringes 5ml, 2-part       | 2005     | 100 pcs / bx<br>20 bx / cs  |
| Disposable Needles<br>for Vaccine Dilution       | ONE-CARE <sup>®</sup> Needles 21G x 1.5"(38mm)             | 21150-CM | 100 pcs / bx<br>100 bx / cs |
| Disinfection                                     | ONE-CARE <sup>®</sup> Alcohol Prep Pads                    | 1021     | 100 pcs / bx<br>100 bx / cs |

WWW.ONE-CARE.COM



MediVena

 Genlake Parkway, Ste 700
 Atlanta, GA 30328, USA
 +1 404 514 2586
 fin@emedivena.com
 www.medivena.com
 CUSTOMER SERVICE TOLL
 FREE NUMBER: 888 270 645

#### Distributed in EU by:

MediVena Sp. z o.o. Al. Jana Pawla II 27 00-867 Warsaw, Poland
 +48 661 920 142
 info@medivena.com www.medivena.com

