

ONE-CARE®

SYRINGES WITH SAFETY NEEDLES



ONE-CARE® Luer-Lock syringes with safety needles help to reduce risk of accidental needle-stick, injuries and transmission of blood borne pathogens for safety of medical personnel and patient.

ONE-CARE® product portfolio includes also 1ml Luer-Lock syringes, 21-23G hypodermic needles and 18G blunt fill needles for vials.

ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 23G x 1" (25mm)


REF 2110 [1 ml] 23G x 1" (25mm)




ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 25G x 1" (25mm)

REF 2112 [1 ml] 25G x 1" (25mm)


Protective shield activation methods:



THUMB ACTIVATION




FINGER ACTIVATION



SURFACE ACTIVATION

Key features:

- Easy one-hand technique of activating safety shield.
- Audible click indicates safety mechanism activation.
- Locking mechanism ensures needle is secure and completely locked in the safety shield.



ONE-CARE®
LUER-LOCK
"ZERO-DEAD SPACE"

MediVena's **ONE-CARE® Syringes and Needles** have "low-dead space" designs to reduce waste of vaccination and to allow for extractions as many doses as possible to maximize the vaccination opportunities for population. In particular, **ONE-CARE®** injection syringes of 1 ml have elongated piston tip design to leave low-to-zero "dead space" within syringe.

ONE-CARE® Safety Needles

Safety Needles 23Gx1" (25mm)

REF 2111



Safety Needles 25Gx1" (25mm)

REF 2113



Key features:

- Reliable safety mechanism integrated with the needle
- Convenient, one-hand use
- Compatible with all Luer-Lock and Luer-Slip Syringes as per ISO 80369
- Bevel oriented needle for easy injection
- Compliant with performance requirements of ISO 7864
- Color coding according to ISO 6009

Before Use




After Use



ONE-CARE® LUER-LOCK DISPOSABLE SYRINGES, 3-PART

REF 2002

1 ml



Luer-Lock 1ml syringes for secure connection with hypodermic needles

Silicone lubrication

Ensures smooth plunger movement

High-transparency barrel

Makes it easy to see syringe content

Safe plunger backstop

Prevents accidental withdrawal of plunger

Large and stable plunger grip

Ensures stability and comfort during injection

ONE-CARE® DISPOSABLE NEEDLES

ONE-CARE® DISPOSABLE NEEDLES 23G X 1" (25mm)

REF 23100-CM

ONE-CARE® DISPOSABLE NEEDLES 23G X 1.5" (38mm)

REF 23150-CM



ONE-CARE® hypodermic needles for vaccine administration have precise tri-bevel design, are electro-polished and lubricated for gentle and virtually pain-free injections.

ONE-CARE® DISPOSABLE NEEDLES 21G X 1.5" (38mm)

REF 21150-CM



ONE-CARE® Canula Needles 21G x 1.5" (38mm) are recommended for vaccine preparation as may be required by vaccine manufacturer's instructions for use.

ONE-CARE® Blunt Fill Needles 18G x 1.5" (38mm) without filter

REF 18150-KD



ONE-CARE® Blunt Fill Needles 18G x 1.5" (38mm) with filter

REF 18150F-KD

ONE-CARE® ALCOHOL PREP PADS

REF 1020




Soft, absorbent, nonwoven pads are saturated with 70% Isopropyl Alcohol which exceeds CDC's recommended 60% alcohol content for optimal antiseptic action.

Intended for Individual, Professional and Hospital Use in preparation of skin prior to an injection or blood collection, cleaning of minor cuts and scrapes, or cleaning of medical equipment and supplies

Compact Box of 100 pads for optimized storage and shipping cost

MEDIUM

▶ **70% Isopropyl Alcohol**
For External Antiseptic Use Only



ORDERING INFORMATION

PRODUCT/APPLICATION	PRODUCT DESCRIPTION	REF	BOX / CASE
Syringes with Safety Needles for Injections	ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 23G x 1" (25mm)	2110	100 pcs / bx 12 bx / cs
	ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 25G x 1" (25mm)	2112	
Safety Needles for Injections	ONE-CARE® Safety Needles 23G x 1" (25mm)	2111	100 pcs / bx 20 bx / cs
	ONE-CARE® Safety Needles 25G x 1" (25mm)	2113	
Syringes	ONE-CARE® Luer-Lock Disposable Syringes, 3-Part	2002	100 pcs / bx 24 bx / cs
Hypodermic Needles	ONE-CARE® Hypodermic Needles 23G x 1" (25mm)	23100-CM	100 pcs / bx 100 bx / cs
	ONE-CARE® Hypodermic Needles 23G x 1.5" (38mm)	23150-CM	
	ONE-CARE® Hypodermic Needles 21G x 1.5" (38mm)	21150-CM	
Blunt Fill Needles for Vials	ONE-CARE® Blunt Fill Needles 18G x 1.5" (38mm) without filter	18150-KD	100 pcs / bx 100 bx / cs
	ONE-CARE® Blunt Fill Needles 18G x 1.5" (38mm) with filter 1.5" (38mm)	18150F-KD	
Disinfection	ONE-CARE® Alcohol Prep Pads	1020	100 pcs / bx 100 bx / cs

WWW.ONE-CARE.COM

Distributed in the US by:

Distributed in EU by:

