

# ONE-CARE<sup>®</sup>

## BLOOD COLLECTION TUBE HOLDERS



Access Holder

Blood Transfer Device



Standard Tube Holder



Safety Holder

### Features

- Size: 0.9 x 19 mm
- Transparent material for better visibility of collection tubes filling up with blood
- Universal design for compatibility with blood collection needles and winged sets
- Tight fit connectors to prevent leakage
- DEHP-free & Latex-free materials

### Sterilization / Expiry

- ETO
- Shelf life: 5 years from mfg date

### Regulatory Compliance

- US FDA, CE, ISO 13485

### Materials

- Medical grade PP
- SUS304
- ABS
- All Materials are latex & DEHP free

## ORDERING INFORMATION

CODE	DESCRIPTION - NEEDLE (DIAMETER x LENGTH)	COLOR	PACKAGING
5015	ONE-CARE <sup>®</sup> Tube Holder	Transparent	100 pcs / box 20 bx (2000 pcs) / case
5016	ONE-CARE <sup>®</sup> Safety Tube Holder	Transparent/Pink	100 pcs / box 10 bx (1000 pcs) / case
5017	ONE-CARE <sup>®</sup> Blood Transfer Device with Female Luer-Lock and Back-End Needle	Transparent/Pink	200 pcs / box 4 bx (800 pcs) / case
5018	ONE-CARE <sup>®</sup> Access Holder with Male Luer Adapter and Back-End Needle	Transparent/Blue	200 pcs / box 4 bx (800 pcs) / case

[WWW.ONE-CARE.COM](http://WWW.ONE-CARE.COM)

**MediVena<sup>®</sup>**  
inspired by innovation

[www.medivena.com](http://www.medivena.com)

**AMERICAS:**  
MediVena  
1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA  
+1 404 514 2586 © info@medivena.com  
**CUSTOMER SERVICE TOLL FREE: 888 270 6457**

**EUROPE:**  
MediVena Sp. z o.o.  
Al. Jana Pawla II 27, 00-867 Warsaw, POLAND  
+48 695 995 700

