

# Lancets and Lancing device



ONE-CARE™ Lancing Devices and Lancets have unique features for easy and seamless blood sampling to enhance patient comfort. ONE-CARE™ Lancing Device has multiple penetration depth settings to adjust to the patient skin type and blood flow requirement. The device has a color-coded loading indicator to easily distinguish when the device is ready to use. After puncture, you can simply slide the lancet ejector forward to safely remove the lancet from the device. ONE-CARE™ Lancets are universal fit and can be used not only with ONE-CARE™ device but most other devices available in the market nowadays. Our needle offering include ultra-thin 32G lancets (~0.2mm diameter almost as thin as a human hair) for comfortable blood collection of tiny drop for modern glucometers.

Color	Specification	Package
Light Blue	23G	100pcs/box 10,000pcs/case
Light Green	26G	
Green	28G	
Purple	30G	
Light Cyan	31G	
Dark Green	32G	

- Universal Lancing Devices
- Mini version: compact size of 9 cm length
- Loading indicator
- Adjustable tip - 7 depth penetration settings for best results.
- Reduced pain - thanks to the linear tracking system
- Quiet and discrete activation
- Safe ejection system - used lancet is easily removed by ejector with just one hand
- 23G-32G Universal fit Lancets

## United States of America

| Atlanta GA (Headquarter)

Address: 1 Glenlake Parkway, Ste 700

Atlanta, GA 30328, USA

Email: info@medivena.com

Phone: (+1) 404 514 2586

## Europe

| Poland (MediVena Sp. Zo.o)

Address: Al. Jana Pawla II 27, 00-867 Warsaw, Poland

Email: info@medivena.com | Phone: (+1) 404 514 2586

| Portugal | Email: info@medivena.com