



Medical Face Mask

Type IIR - EN 14683

Medical Face Mask

Type IIR - EN 14683

(BFE)
 $\geq 98\%$
 Bacterial Filtration Efficiency



Medical Face Mask with Ear Loops

EN 14683 Type IIR

Single-Use

Fluid Protection

Soft & Easy to breathe through

Other features:


- Size 17 x 9,5 cm
- Non-Sterile
- Not made with natural rubber latex

The table below indicates the performance requirements defined by **EN 14683-2019**. This product meets European Standard **EN 14683** for medical face masks **Type IIR**. The product is also covered by **Declaration of Conformity** with European Directive, Document No. **KB/CE-011-03**. In addition, the manufacturer's Quality Management System is certified for compliance with **EN ISO13485:2016** by TÜV SÜD Product Service GmbH, Germany, Certificate No. Q6 050970 0016 Rev. 01.

Conformity with European Norms, Medical Device Directives and International Standards ISO

ORDERING INFORMATION

Test	Type / Typ I *	Type / Typ II	Type / Typ IIR
Bacterial filtration efficiency (BFE) (%)	$\geq 95\%^*$	$\geq 98\%$	$\geq 98\%$
Differential pressure (Pa/cm ²)	$< 40\%^*$	< 40	< 60
Splash resistance pressure (kPa)	Not required*	Not required*	$\geq 16,0$
Microbial cleanliness (cfu/g)	$\leq 30\%^*$	≤ 30	≤ 30

REF	 KB/IIR
Product Description	Medical Face Mask with Ear Loops, Type IIR, Box of 50
Box Quantity	Box of 50 pcs (5 Bags of 10 pcs)
Case Quantity	2,000 pcs (40 boxes of 50 pcs)
Dimensions	55,5 * 44,5 * 42cm
Weight	11 kg

*Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.



CUSTOMER SERVICE TOLLFREE NUMBER

888 270 6457

MediVena
 inspired by innovation

Distributed in the US by:
MediVena

1 Glenlake Parkway, Ste 700
 Atlanta, GA 30328, USA
 +1 404 514 2586
 info@medivena.com

Distributed in European Union by:
Medivena Sp. z o.o.

Aleja Jana Pawła II 27
 00-867 Warszawa, POLAND
 +48 661 920 142
 info@medivena.com

